



more consistent yields


high-performance
disposables for pcr

FROM BRIDGE BIOSCIENCE CORPORATION

- reduced wall thickness
- better thermal response
- faster cycling
- more uniform yields

The extra efficiency and uniformity of UTW plates and tubes can improve pcr yields and reduce "drop outs" caused by inconsistent thermal conditions. These results are particularly noticable during "fast" pcr protocols.

UTW = more uniform reactions



Comparison of UTW and competitor tubes run side by side during identical 2-step high-speed pcr reactions. Samples from all reactions were run on an agarose gel to compare band strength of the 400bp reaction product.

UTW & Ultra-ThinWall are trademarks of BioInnovations Oy. US and international patents pending. 02/2007

BRIDGE

BIOSCIENCE CORPORATION

www.bridgebio.com

200 International Dr Suite 125
Portsmouth NH 03801
603-766-5990 info@bridgebio.com